

## Standard Operating Procedure

<b>Title</b>	<b>Viability protocol by using a Cell Counter</b>
Subtitle	
NANoREG Work package/task:	WP5 / Task 5.06
Owner and co-owner(s)	Ricard Marcos, UAB
Date finalised	June 2017
Document name	NANoREG D5.07 SOP 05 Viability protocol by using a Cell Counter
Key words:	

<b>Version</b>	<b>Date</b>	<b>Reason of change</b>
1	2017/07/10	Harmonisation according NANoREG template for SOPs
2		
3		
4		

This work is licensed under the Creative Commons Attribution-NonCommercial-ShareAlike 4.0 International License.

To view a copy of this license, visit <http://creativecommons.org/licenses/by-nc-sa/4.0/> or send a letter to Creative Commons, PO Box 1866, Mountain View, CA 94042, USA.

## 1 Materials

- Desired cells
- PBS
- PBS + 0.1% FBS
- Trypsin (in case of adherent cultures)
- Cell culture medium
- Beckman-Coulter Cell Counter

## 2 Procedure

- 1) In the case of adherent cultures, trypsinize and resuspend the cells in a known volume of PBS + 0.1% FBS. The volume of trypsin, its concentration and the trypsinization time depends on the cell line. Cells that grow in suspension can be counted directly.
- 2) Prepare a cuvette with 9.9 mL of PBS and add 100  $\mu$ L of your cell culture (1:100 dilution).
- 3) Set the Coulter Counter:
  - a. Set the dilution factor to 1 E + 02.
  - b. Set the approximate size of the cells, entering the lower limit (TI) and higher (You).
- 4) Shake the sample without making bubbles. Place the electrode tube inside the cuvette.
- 5) Press start. When the machine has finished, press "output" and manually adjust your cells' size to refine the count. This gives you the amount of cell/mL. Multiply this value for the volume you have, to get the total number of cells you have in your cell culture.