

## Standard Operating Procedure

<b>Title</b>	<b>Confocal Raman Microspectroscopy (CRM)</b>
Subtitle	
NANoREG Work package/task:	WP5 / Task 5.06
Owner and co-owner(s)	Irina Estrela-Lopis, ULEI Carolin Merker, ULEI
Date finalised	June 2017
Document name	NANoREG D5.07 SOP 18 Confocal Raman Microspectroscopy
Key words:	

<b>Version</b>	<b>Date</b>	<b>Reason of change</b>
1	2017/07/10	Harmonisation according NANoREG template for SOPs
2		
3		
4		

This work is licensed under the Creative Commons Attribution-NonCommercial-ShareAlike 4.0 International License.

To view a copy of this license, visit <http://creativecommons.org/licenses/by-nc-sa/4.0/> or send a letter to Creative Commons, PO Box 1866, Mountain View, CA 94042, USA.

## Table of content

OBJECTIVE.....	2
PRINCIPLE OF THE METHOD.....	2
EQUIPMENT & CHEMICALS.....	2
SETTINGS FOR CRM MEASUREMENTS.....	3
SAMPLE PREPARATION.....	3
DATA ANALYSIS.....	3

### Objective

Confocal Raman Microspectroscopy (CRM) as a label-free and non-invasive imaging technique provides unique information about translocation and fate of nanoparticles at subcellular level. CRM allows the simultaneous visualization of NMs and their biological environment. CRM has a strong potential as a label-free, nondestructive technique for time-course imaging of individual cells and tracking of cell metabolism.

### Principle of the method

The combination of confocal microscopy with Raman spectroscopy provides very high content of molecular-spectroscopic information in every pixel of image at a singular cell level. CRM offers a label-free way to complement histopathology, diagnose and characterize culture cells and tissue sections. CRM provides 3D chemical composition images with a resolution of about 250 nm. CRM reveals not only the 3D NM distribution but also their co-localization within cell compartments.

### Equipment & chemicals

1. Confocal Raman Microspectroscopy (Alpha300 R, WITec GmbH, Germany)
2. 532, 630, and 785 nm laser source
3. 63x water immersion objective (W Plan-Apochromat 63x/1, Zeiss, Germany)
4. Heating stage (WITec GmbH, Germany)
5. CaF<sub>2</sub> substrate: 25mm x 25mm (CRYTAL GmbH, Germany), quartz slices: 25mm diameter (ZELL QUARZGLAS und Technische Keramik Technologie GmbH, Germany), silica wafers
6. Superfibronectin from human plasma (S5171, Sigma, St. Luis, MO)
7. Dulbecco`s Phosphate Buffered Saline (PBS w/o Ca/Mg) (Biowest, France)
8. Formaldehyde
9. Xylol

## Settings for CRM measurements

- Laser: 532nm
- Grating: 600 g mm<sup>-1</sup>
- Laser power: NM212, NM103, NM104, NM110 & lung tissue: 34mW  
NM300-K: 16mW
- Integration time: NM212, NM103, NM104, NM110 & lung tissue: 0.07s  
NM300-K: 0.14s

## Sample preparation

### 1. Culture cells

- Cell cultivation on CaF<sub>2</sub>, quartz or silica wafer (NP dependent)
- Coating of substrate with superfibronectin (5µg/ml in PBS): 2h at 37°C
- Wash 3-times with PBS, cultivate cells in culture medium for 24h (cell density: 75000 cells/ml)
- NP exposure:
  1. *during the exposure*: live cell imaging of NPs
  2. *after the exposure*: remove cell medium with NPs, wash 3-times with PBS, fix cells 30min at room temperature with 3.7% Formaldehyde, wash 3-times with PBS, store cells at 4°C and measure in PBS

### 2. Lung slides

- Paraffin-embedded tissues on quartz slides: removal of paraffin by means of xylol (3-times, 10min)
- Measure in Millipore-water

## Data analysis

Acquired spectra were processed using the Project FOUR PLUS software (WITec GmbH, Germany).