

"SPOLIGOTYPING"

a PCR-based method to simultaneously detect and type
Mycobacterium tuberculosis complex bacteria

**Kristin Kremer, Annelies Bunschoten, Leo Schouls,
Dick van Soolingen, and Jan van Embden**

**National Institute of Public Health and the Environment,
Bilthoven, The Netherlands**

Version November 2002

CONTENTS

	page
1. General description	
3	
1.1. Principle of the spoligotyping method	
3	
1.2. Practical use	6
2. Methods	7
2.1. Preparation of the membrane containing the spacer-oligonucleotides	7
2.2. <i>In vitro</i> amplification of spacer DNA by PCR	8
2.3. Hybridization with PCR product and detection	9
2.4. Reuse of the membrane	10
3. Solutions	11
4. Supplies	12
5. Troubleshooting	13
6. References	14
7. Supplements	16
1. Concentration of the aminolink-oligo's	16
2. Template DNA isolation methods	17
3. Spacer sequences	19

1. GENERAL DESCRIPTION

1.1 Principle of the spoligotyping method

The typing method described in this protocol is based on DNA polymorphism present at one particular chromosomal locus, the "Direct Repeat" (DR) region, which is uniquely present in *Mycobacterium tuberculosis* complex bacteria. This locus was first described by Hermans et al.⁵ who sequenced this region in *Mycobacterium bovis* BCG, the strain used worldwide to vaccinate against tuberculosis. The DR region in *M. bovis* BCG consists of directly repeated sequences of 36 base pairs, which are interspersed by non-repetitive DNA spacers, each 35 to 41 base pairs in length. The number of copies of the DR sequence in *M. bovis* BCG was determined to be 49. In other *M. tuberculosis* complex strains the number of DR elements was found to vary significantly. The vast majority of the *Mycobacterium tuberculosis* strains contain one or more IS6110 elements in the DR region (Fig.1).




Figure 1. Structure of the DR locus in the genome of *M. tuberculosis* H37Rv and *M. bovis* BCG P3. The rectangles depict the 36 bp Direct Repeat (DR).

In contrast to the DRs, the spacers are usually present only once in the DR region, but occasionally some are found twice, either separated by one or by several DR's and other spacers. One DR and its neighbouring non-repetitive spacer is termed "Direct Variant Repeat" (DVR).

When the DR regions of several strains were compared, it was observed that the order of the spacers is about the same in all strains, but deletions and/or insertions of spacers and DR's occur (Fig. 2)². The mechanism by which spacers and copies of DR are generated, is unknown. With the method described here, the presence or absence in the DR region of 43 spacers of known sequence can be detected by hybridization of PCR-amplified spacer DNA to a set of immobilized oligonucleotides, representing each of the unique spacer DNA sequences. This method will be referred to as *spoligotyping* (from *spacer oligotyping*)⁷.




Figure 2. Schematic presentation of the polymorphism in DR regions of different *M. tuberculosis* complex strains. Blocks of DVR are missing in one strain when compared to another. The spacer order remains about the same.

By spoligotyping one can detect the presence or absence of spacers of known sequence. The first step in the method is to amplify the DR region of a given strain by PCR. The primers used are based on the sequence of the DR, and allow the amplification of the spacer(s) between the DR targets (Fig. 3). The obtained PCR products differ in length because of two reasons. First, the product contains several spacers and the DR's in between if the primers anneal to DR's not next to each other. Second, the product itself can act as a primer, and become elongated with one or more DVRs. Therefore, the PCR product provides no reliable information about spacer order or total length of the DR region. A biotin labelled reverse primer is used, so that all the reverse strands synthesized are biotin labelled.




Figure 3. Principle of the *in vitro* amplification of DNA within the DR region of *M. tuberculosis* complex bacteria. The use of the 2 primers, a and b, for *in vitro* amplification, will lead to the amplification of any spacer or a stretch of neighbouring spacers and DR's.

Oligonucleotides derived from the known spacers in the DR cluster are covalently linked to an activated Biodyne C membrane in parallel lines. PCR products are hybridized perpendicular to the oligo lines. After hybridization the membrane is incubated in streptavidin peroxidase, which binds to the biotin label on the PCR products. Detection of hybridization signals is done by the enhanced chemiluminescence (ECL) detection system. The peroxidase present on the streptavidine catalyzes a reaction resulting in the emission of light, which can be detected by autoradiography of the membrane. We refer to this type of blot as *reversed line blot* (Fig. 4).

Spoligotyping




Figure 4. Overview of the spoligotyping method.

An example of a result of the spoligotyping method used to analyze a variety of clinical isolates is shown in Fig. 5.




Figure 5. Spoligotyping result of *M. tuberculosis* H37Rv, *M. bovis* BCG P3 and 38 different clinical isolates. A membrane with 43 spacer oligonucleotides was used (vertical lines). The spacer oligonucleotides were derived from the spacers of *M. bovis* BCG P3, *M. tuberculosis* H37Rv.

1.2 Practical use

Spoligotyping may offer an alternative for typing Southern blotting when rapid results are required. The method is in particular useful to simultaneously detect and type *M. tuberculosis* complex bacteria in clinical samples (suspected nosocomial infections, outbreaks in prisons, etc.). The level of differentiation by spoligotyping is less compared to IS6110 fingerprinting for strains having five or more IS6110 copies, but higher for strains with less than five copies. Thus spoligotyping is a preferred method to type *M. bovis* strains, which usually contain only one or two IS6110 copies. Note that *M. bovis* can be recognized by the absence of reactivity with spacers 39-43 (Fig. 5).

2. METHODS

2.1 Preparation of the membrane containing the spacer-oligonucleotides

Purpose

Preparation of a membrane containing the spacer-oligonucleotides. The spacer-oligonucleotide sequences are derived from DNA sequences of the DR region in the strains *M. tuberculosis* H37Rv and *M. bovis* BCG P3. All spacer-oligonucleotides are synthesized with a 5' terminal aminogroup, by which they can be covalently linked to an activated negatively charged Biodyne C membrane.

Procedure

- 1) Dilute the spacer oligonucleotides to the optimized concentrations (see supplement 1) in 150 μ l 500 mM NaHCO₃, pH 8.4.^{a)}
- 2) Activate the Biodyne C membrane by 10 min incubation in 10 ml freshly prepared 16% (w/v) 1-ethyl-3-(3-dimethylaminopropyl) carbodiimide (EDAC) in demineralized water, in a rolling bottle at room temperature.
- 3) Place the membrane in a plastic container and shake with demineralised water for 2 min and place it on a support cushion in a clean miniblotted system. Turn the screws hand-tight.
- 4) Fill the slots of the miniblotted with 150 μ l of the diluted oligonucleotide solutions but do not use the first and the last slot to apply oligo's.
- 5) Use the first and the last slot to mark the edges of the membrane by adding drawing pen ink, diluted 1:100 in 2xSSPE buffer.
- 6) When all samples are added, incubate for 2 min at room temperature.
- 7) Remove the oligonucleotide solutions by aspiration in the same order as in which they were applied to the slots.
- 8) Remove the membrane from the miniblotted using forceps and incubate it in 250 ml freshly prepared 100 mM NaOH for 10 min (maximum) in a plastic container, while shaking, to inactivate the membrane.^{b)}
- 9) Rinse the membrane with demineralised water.
- 10) Wash the membrane in a plastic container by gentle shaking in 250 ml 2 x SSPE/0.1% SDS for 5 min at 60 C.
- 11) Wash the membrane in a plastic container by gentle shaking in 100 ml 20 mM EDTA pH 8 for 15 min at room temperature.
- 12) Store the membrane at 4 C until use (sealed in plastic or wrapped in Saran-wrap, to avoid dehydration of the membrane).

Remarks

- a) The pH of the NaHCO₃ should be **exactly** 8.4.
- b) Inactivation of the membrane by NaOH longer than 10 min will result in weak hybridization signals.

2.2 *In vitro* amplification of spacer DNA by PCR

Principle

Amplification of the spacers is accomplished by using the primers DRa and DRb, which enable to amplify the whole DR region. Only a very small amount of template DNA is required. Typically the PCR is performed on 10 ng purified chromosomal Mycobacterial DNA but, with minor adaptations, DNA extracts from clinical samples or lysed bacteria (from freezer or Löwenstein) can also serve as template (see supplement 2). The PCR products are labeled with biotin, because one of the primers is biotinylated. The primers for the PCR are based on the DR sequence:

DRa : 5' -GGT TTT GGG TCT GAC GAC-3' , biotinylated at 5' end.

DRb : 5' -CCG AGA GGG GAC GGA AAC-3'

Procedure

- 1) Dilute the DNA samples to the required concentration. Always include chromosomal DNA of *M. tuberculosis* strain H37Rv and *M. bovis* BCG as positive controls. Use water as a negative control.
- 2) Prepare the reactionmixture: ^{a)}
 - x μ l template DNA(10 ng)
 - 4 μ l primer DRa (20 pmol) ^{b)}
 - 4 μ l primer DRb (20 pmol) ^{c)}
 - 4 μ l dNTP-mixture (2.5 mM each dNTP, final conc. 0.2 mM each dNTP)
 - 5 μ l 10x concentrated Super *Tth* buffer
 - 0.1 μ l Super *Tth* polymerase (5 units/ μ l)
 - 33-x μ l MQ water (to a final volume of 50 μ l)
- 3) Add one drop of mineral oil to the tubes to prevent evaporation of the PCR-mix during the amplification.
- 4) Place the tubes in a PCR-apparatus for amplification, and perform the following temperature cycling: ^{d,e,f)}
 - 3 min 96 C
 - 1 min 96 C
 - 1 min 55 C 20x
 - 30 sec 72 C
 - 5 min 72 C

Remarks

- a) Preparation of the mixture has to take place in a laboratory free of mycobacterial PCR products containing the DR sequences.
- b) Primer DRa is biotinylated and should be stored at +4°C. Repeated freeze-thawing of the biotinylated primer results in weaker Spoligopatterns.
- c) Primer DRb should be stored in small aliquots at -20 C.
- d) For amplification of the DR-cluster from extracts of clinical samples, the number of cycles can be increased to 40.
- e) For amplification of the DR-cluster from heat-killed cells, the number of cycles can be increased to 30 (see also supplement 2B, Preparation of lysates from colonies).
- f) PCR products can be used immediately, but can also be stored at -20°C to be used later.

2.3 Hybridization with PCR product and detection

Purpose

Hybridization of the biotin-labeled PCR products to the immobilized spacer-oligos that represent spacers of known sequence. The presence of spacers is visualized on film as black squares after incubation with streptavidin-peroxidase and ECL-detection.

Note

All incubations should take place in a plastic container under gentle shaking, unless otherwise stated.

Thoroughly clean the miniblottedter with soap and a dedicated brush before use.

Never touch the membrane with gloves, the powder causes background. Use forceps.

The quality of the SDS is of critical importance. It should be fresh, do not store it for longer than one week. We have good experience with SDS from BDH Laboratory Supplies.

Procedure

- 1) All buffers should be prewarmed before use. Prepare the following buffers from concentrated stocks, using demineralized water for dilution (quantities for one membrane):
 - 250 ml 2xSSPE/0.1% SDS, 60°C,
 - 500 ml 2xSSPE/0.5% SDS, 60°C,
 - 500 ml 2xSSPE/0.5% SDS, 42°C.
 - 500 ml 2xSSPE, room temperature.
- 2) Add 20 µl of the PCR products to 150 µl 2xSSPE/0.1% SDS.
- 3) Heat-denature the diluted PCR product for 10 min at 99 C and cool on ice immediately.
- 4) Wash the membrane for 5 min at 60 C in 250 ml 2xSSPE/0.1% SDS.
- 5) Place the membrane and a support cushion into the miniblottedter, in such a way that the slots are perpendicular to the line pattern of the applied oligonucleotides. ^{a)}
- 6) Remove residual fluid from the slots of the miniblottedter by aspiration.
- 7) Fill the slots with the diluted PCR product (avoid air bubbles!) and hybridize for 60 min at 60°C on a horizontal surface (no shaking!). Avoid contamination of neighbouring slots. ^{b)}
- 8) Remove the samples from the miniblottedter by aspiration and take the membrane from the miniblottedter using forceps.
- 9) Wash the membrane twice in 250 ml 2xSSPE/0.5% SDS for 10 min at 60 C .
- 10) Place the membrane in a rolling bottle and allow it to cool down to prevent inactivation of the peroxidase in the next step.
- 11) Add 2.5 µl streptavidin-peroxidase conjugate (500U/ml) to 10 ml of 2xSSPE/0.5% SDS, and incubate the membrane in this solution for 45 to 60 min at 42 C in the rolling bottle.
- 12) Wash the membrane twice in 250 ml of 2xSSPE/0.5% SDS for 10 min at 42 C.
- 13) Rinse the membrane twice with 250 ml of 2xSSPE for 5 min at room temperature.
- 14) For chemiluminiscent detection of hybridizing DNA, incubate the membrane for 1 min in 20 ml ECL detection liquid. ^{c)}

- 15) Cover the membrane with a transparent plastic sheet or Saran-wrap and expose a light sensitive film to the membrane for 20 min. ^{d)}
- 16) If the signal is too weak or too strong the membrane can be used again directly to expose another film for a shorter or longer period.

Remarks

- a) Do not reuse the supportcushions.
- b) If less than 45 samples are applied to the miniblottedter, fill one neighbouring slot with 2XSSPE/0.1%SDS to prevent cross-flow.
- c) Use a dedicated plastic container. Do not use this container for other purposes, since some reagents decrease the intensity of the Spoligopatterns.
- d) If the result is unsatisfactory, you can try to improve this.
Black spots (background) possibly occur due to contamination during filter handling (e. g. touched with fingers). Start again from step 8.
Blank areas in the spoligopatterns possibly indicate that the membrane was not completely soaked with ECL detection liquid. Start again from step 13.

2.4 Reuse of the membrane

Purpose

The hybridized PCR product is dissociated from the membrane in order to reuse the membrane for the next hybridization. A membrane can be reused for about 15 times.

Procedure

- 1) Wash the membrane twice by incubation in 1% SDS at 80 C for 30 min.
- 2) Wash the membrane in 20 mM EDTA pH 8, for 15 min at room temperature.
- 3) Store the membrane at 4 C until use (sealed in plastic or wrapped in Saran-wrap, to avoid dehydration of the membrane).

3. SOLUTIONS

500 mM NaHCO₃, pH8.4

10.5 g NaHCO₃ / 250 ml demineralised water.

16% (w/v) EDAC

1.6 g EDAC / 10 ml demineralised water.

100 mM NaOH

0.8 g NaOH / 200 ml demineralised water.

0.5 M EDTA, pH 8.0

186.12 g EDTA / L demineralised water.

20 mM EDTA, pH 8.0

Dilute 0.5 M EDTA 25 times.

20xSSPE

0.2 M Na₂HPO₄*2H₂O 35.6 g/l

3.6 M NaCl 210.24 g/l

20 mM EDTA 7.4 g/l

The pH should be 7.4. Autoclave.

Store at roomtemperature for no longer than one year.

2xSSPE

Dilute 20xSSPE ten times with demineralised water.

10% SDS

10 g SDS / 100 ml demineralised water.

2xSSPE/0.1%SDS

Add 100 ml 20xSSPE and 10 ml 10% SDS to 890 ml demineralised water.

2xSSPE/0.5%SDS

Add 100 ml 20xSSPE and 50 ml 10% SDS to 850 ml demineralised water.

4. SUPPLIES

Miniblotter MN45	± 1,150 USD
Order number: MN45 Immunitics, 380 Green street, Cambridge, Mass. 02139, USA Tel. 617-492-5416, Fax. 617-868-7879, Toll-free: 800-227-4765	
Foam cushions	± 200 USD
Order number: PC200 Immunitics, 380 Green street, Cambridge, Mass. 02139, USA Tel. 617-492-5416, Fax. 617-868-7879, Toll-free: 800-227-4765	
Biodyne C membrane	± 150 USD
Order number: BNBCH5R Pall Biosupport	
EDAC (100 g)	± 525 USD
Order number: E.7750 SIGMA	
20 x SSPE, 4 ltr.	± 125 USD
Order number: 15591-035 Gibco BRL Life Technologies Inc.	
SDS specially pure (500 g)	± 580 USD
Order number: 44244 4H BDH Laboratory Supplies	
Super Tth buffer (5 ml)	± 25 USD
Order number: TPRB SphaeroQ	
Super Tth DNA polymerase (5000 U)	± 1,750 USD
Order number: TPO3C SphaeroQ	
Streptavidin-POD-conjugate	± 225 USD
Order number:1089153 Boehringer	
ECL detection liquid	± 1,000 USD
Order number: RPN2105 Amersham International	
Hyperfilm ECL	± 40 USD
Order number: RPN 2103 Amerham International	

5. TROUBLESHOOTING

- 1) No hybridization signal detected: analyze 5 µl of the PCR product on a 2% agarose gel. A ladder pattern should be visible. If a ladder pattern is visible, check the labelling of the PCR product by spotting it onto a membrane, followed by incubation with streptavidin peroxidase.
- 2) High background (stripes): thoroughly clean the miniblotter using a dedicated brush, and soak the apparatus, preferably overnight, in a soap solution, e.g. Extran (Merck).
- 3) High background (spots): strip the membrane again, and test it with PCR products of the control strains. If stripping does not lead to a lower background, the membrane should not be used anymore.
- 4) In case of other problems, contact us, preferably by e-mail:

Kristin Kremer

E: Kristin.Kremer@rivm.nl

T: +31-30-2742720

Dick van Soolingen

E; d.van.soolingen@rivm.nl

T: +31-30-2742363

Mycobacteria Reference Unit (pb 22),
Diagnostic Laboratory for Infectious Diseases and Perinatal Screening,
National Institute of Public Health and the Environment,
P.O. Box 1, 3720 BA, Bilthoven, The Netherlands.

Fax: +31-30-2744418

6. REFERENCES

1. **Aranaz, A., E. Liebana, A. Mateos, L. Dominquez, D. Vidal, M. Domingo, O. Gonzalez, E. F. Rodriguez-Ferri, A. E. Bunschoten, J. D. A. van Embden, and D. Cousins.** 1996. Spoligotyping of *Mycobacterium bovis* strains from cattle and other animals: a tool for epidemiology of tuberculosis. *J. Clin. Microbiol.* **34**:2734-2740.
2. **Groenen, P. M. A., A. E. Bunschoten, D. van Soolingen, and J. D. A. van Embden.** 1993. Nature of DNA polymorphism in the direct repeat cluster of *Mycobacterium tuberculosis*; Application for strain differentiation by a novel method. *Mol. Microbiol.* **10**:1057-1065.
3. **Goguet de la Salmoniere, Y. O., H., M. Li, G. Torrea, A. Bunschoten, J. van Embden, and B. Gicquel.** 1997. Spoligotyping: evaluation of this method by studying the transmission of *Mycobacterium tuberculosis*. Combination with inter-*IS6110* PCR for typing *M. tuberculosis* complex strains. *J. Clin. Microbiol.* **35**:2210-2214.
4. **Goyal, M., N. A. Saunders, J. D. A. van Embden, D. B. Young, and R. J. Shaw.** 1997. Differentiation of *Mycobacterium tuberculosis* isolates by Spoligotyping and *IS6110* restriction fragment length polymorphism. *J. Clin. Microbiol.* **35**:647-651.
5. **Hermans, P. W. M., D. van Soolingen, E. M. Bik, P. E. W. de Haas, J. W. Dale, and J. D. A. van Embden.** 1991. The insertion element *IS987* from *M. bovis* BCG is located in a hot spot integration region for insertion elements in *M. tuberculosis* complex strains. *Infect. Immun.* **59**:2695-2705.
6. **Heyderman, R. S., M. Goyal, P. Roberts, S. Ushewokunze, S. Zizhou, B. G. Marshall, R. Makombe, J. D. A. van Embden, P. R. Mason, and R. J. Shaw.** 1998. The epidemiology of sputum-positive pulmonary tuberculosis in Harare, Zimbabwe - analysis by spoligotyping. *Thorax*: **53**:346-350.
7. **Kamerbeek, J., L. M. Schouls, A. Kolk, M. van Agterveld, D. van Soolingen, S. Kuijper, A. E. Bunschoten, H. Molhuizen, R. Shaw, M. Goyal, and J. D. A. van Embden.** 1997. Simultaneous detection and strain differentiation of *Mycobacterium tuberculosis* for diagnosis and epidemiology. *J. Clin. Microbiol.* **35**:907-914.
8. **Kaufhold, A., A. Podbielski, G. Baumgarten, M. Blokpoel, J. Top, and L. Schouls.** 1994. Rapid typing of group A streptococci by the use of DNA amplification and nonradioactive allele specific oligonucleotide probes. *FEMS Microbiol. Lett.* **119**:19-26.

9. **Mangiapan, G., M. Vokurka, L. Schouls, J. Candranel, D. Lecossier, J. van Embden, and A. Hance.** 1996. Sequence capture - PCR improves the detection of mycobacterial DNA in clinical specimens. *J. Clin. Microb.* **34**:1209-1215.
10. **Samper, S., C. Martin, A. Pinedo, A. Rivero, J. Blázquez, F. Baquero, D. van Soolingen, and J. van Embden.** 1997. Detection of international transmission of multidrug-resistant tuberculosis by using DNA fingerprint databases. *AIDS*: **11**:1237-1242.
11. **Van Embden, J. D. A., T. van Gorkom, K. Kremer, R. Jansen, B. A. M. van der Zeijst, and L. M. Schouls.** 2000. Genetic variation and evolutionary origin of the direct repeat locus of *Mycobacterium tuberculosis* complex bacteria. *J. Bact.* **182**:2393-2401.
12. **Van Soolingen, D., L. Qian, P. E. W. de Haas, J. T. Douglas, H. Traore, F. Portaels, H. Z. Qing, D. Enkhasaikhan, P. Nymadawa, and J. D. A. van Embden.** 1995. Predominance of a single genotype of *Mycobacterium tuberculosis* in countries of East Asia. *J. Clin. Microb.* **33**:3234-3238.

SUPPLEMENT 1**Concentration of the aminolink-oligo's.**

Oligo number	Oligo name	Concentration (pmol/150 µl)	Oligo number	Oligo name	Concentration (pmol/150 µl)
1	SP-F310A	12.5	23	SP-F510A	50.0
2	SP-F320A	30.0	24	SP-F520A	50.0
3	SP-F330A	12.5	25	SP-F530A	25.0
4	SP-F340A	12.5	26	SP-F540A	12.5
5	SP-F350A	12.5	27	SP-F550A	25.0
6	SP-F360A	12.5	28	SP-F560A	12.5
7	SP-F370A	25.0	29	SP-F570A	12.5
8	SP-F380A	50.0	30	SP-F580A	12.5
9	SP-F390A	12.5	31	SP-F590A	12.5
10	SP-F400A	15.0	32	SP-F600A	25.0
11	SP-F410A	30.0	33	SP-F610A	100
12	SP-F420A	60.0	34	SP-F620A	25.0
13	SP-F430A	12.5	35	SP-F630A	12.5
14	SP-F440A	30.0	36	SP-F640A	12.5
15	SP-F450A	30.0	37	SP-F650A	12.5
16	SP-F455A	12.5	38	SP-F660A	25.0
17	SP-F457A	100	39	SP-F670A	25.0
18	SP-F460A	12.5	40	SP-F680A	12.5
19	SP-F470A	12.5	41	SP-F690A	12.5
20	SP-F480A	12.5	42	SP-F700A	25.0
21	SP-F490A	25.0	43	SP-F710A	50.0
22	SP-F500A	12.5			

SUPPLEMENT 2

Template DNA isolation methods

General remarks:

For handling clinical samples that are to be used in a PCR it is recommended to work in a room specially equipped for this purpose (over-pressure, laminar-flow hood).

Always use negative controls.

A. Preparation of chromosomal DNA

Typically, isolation and purification of chromosomal *Mycobacterium* DNA is done using the CTAB method:

- 1) Transfer at least one loopfull of cells into an Eppendorf tube containing 400 μ l of 1xTE.
- 2) Heat 20 min at 80°C to kill the cells, and cool to room temperature.
- 3) Add 50 μ l 10 mg/ml lysozyme, vortex and incubate at least 1 hr at 37°C.
- 4) Add 75 μ l 10 % SDS/proteinase K solution (5 μ l proteinase K, 10 mg/ml and 70 μ l 10% SDS), vortex shortly and incubate 10 min at 65°C.
- 5) Add 100 μ l 5 M NaCl.
- 6) Add 100 μ l CTAB/NaCl solution (4.1 g NaCl and 10 g CTAB [N-cetyl-N,N,N,-trimethylammoniumbromide] in 100 ml distilled water), which is prewarmed at 65°C. Vortex until the liquid content becomes white ("milky"). Incubate 10 min at 65°C.
- 7) Add 750 μ l of chloroform/isoamyl alcohol (24:1), vortex 10 sec, and centrifuge at room temperature for 5 min, 12,000 *g*.
- 8) Transfer the aqueous supernatant to a fresh microcentrifuge tube.
- 9) Add 450 μ l isopropanol.
- 10) Incubate 10 min on ice.
- 11) Centrifuge 15 min at room temperature.
- 12) Discard the supernatant and wash the pellet with 1 ml of 70% ethanol and centrifuge (approximately 5 min at room temperature).
- 13) Discard the supernatant and dry the pellet.
- 14) Redissolve the pellet in 20 μ l of 1xTE buffer. The DNA can be stored at 4°C until further use.

B. Preparation of lysates from colonies

- 1) Resuspend 2 loops of cells in 250 μ l 1xTE in an Eppendorf tube.
- 2) Kill the cells by incubation at 80 C for 1 hour.
- 3) Centrifuge the tube at 13000 rpm for 2 min, discard the supernatant and resuspend the pellet in 500 μ l of 150 mM NaCl. Repeat this step twice.
- 4) Discard the supernatant and resuspend the pellet in 25 μ l of distilled water or 1x TE.

C. Extraction of total DNA from clinical samples

- 1) Bring the sample (max 1 ml or 1 cm³) aseptically in a 10 ml tube with 2-3 ml of digestion buffer (500 mM Tris HCl, pH 9.0, 20 mM EDTA, 10 mM NaCl, 1%SDS) and incubate overnight at 60 C.
- 2) Vortex the sample for 20 sec, add 0.5 ml phenol to 0.9 ml sample and vortex for 20 sec.
- 3) Centrifuge for 5 min at max speed.
- 4) Transfer the aqueous phase to a fresh tube, containing 0.5 ml phenol, vortex for 20 sec and centrifuge for 5 min at max speed.
- 5) Transfer the aqueous phase (approx. 350 µl) to a fresh tube containing 35 µl 3 M NaAc and 800 µl absolute ethanol, mix and incubate for 20 min at -20 C.
- 6) Centrifuge for 30 min at room temperature, max speed.
- 7) Discard the supernatant and wash the pellet with 500 µl 70% ethanol, centrifuge for 5 min at max speed.
- 8) Discard the supernatant, dry the pellet, resuspend the DNA in 50-200 µl 1x TE and store at -20 C until further use.

D. Isolation of genomic DNA from paraffin-embedded tissues

Sample preparation of DNA from the paraffin-embedded tissues is successful, but remains time consuming. The technique involves two major steps, deparaffinization and protein digestion, each of which involves several centrifugations and washes and requires multiple tube transfers. We use a method which is a one step procedure without protein digestion:

- 1) Add 150 µl of a 5% Chelex suspension to a 14 µm paraffin-embedded tissue section.^{a, b)}
- 2) Vortex thoroughly. The section should be completely covered with the Chelex suspension.
- 3) Heat the mixture 30 min at 100°C. The paraffin then appears floating on the surface of the solution.
- 4) Centrifuge 10 min at 13000 x g.
- 5) Transfer the solution beneath the paraffin, containing the extracted DNA, to a clean microcentrifuge tube.^{c)}
- 6) The PCR is run with two different dilutions of the extracted DNA: 10 µl undiluted DNA and 10 µl DNA of an 1:4 dilution, in 50 µl PCR-mix.

Remarks

- a) To prepare the sections use another scalpel or knife for each sample. The microtome and the knife should be disinfected with 1 N HCl for two minutes after each sample to prevent contamination. Cut a negative control between every sample.
- b) The aim of Chelex®100 treatment is to remove metal ions. Chelex® 100 is stable for at least 2 years when stored sealed in the original container at 22°C. If left in the hydrogen form for more than a few hours, the Chelex has a tendency to lose chelating capacity.
- c) Avoid transferring the Chelex. Chelex will bind Mg²⁺ in the PCR-mix.

SUPPLEMENT 3**Spacer sequences**

The sequences of the spacer-specific oligonucleotides are:

1: ATAGAGGGTCGCCGGT <u>T</u> CTGGATCA ¹⁾	23: AGCATCGCTGATGCGGTCCAGCTCG
2: CCTCATA <u>A</u> ATTGGGCGACAGCTTTTG ²⁾	24: CCGCCTGCTGGGTGAGACGTGCTCG
3: CCGTGCTTCCAGTGATCGCCTTCTA	25: GATCAGCGACCACCGCACCTGTCA
4: ACGTCATACGCCGACCAATCATCAG	26: CTTCAGCACCACCATCATCCGGCGC
5: TTTTCTGACCACTTGTGCGGGATTA	27: GGATTCGTGATCTCTTCCC GCGGAT
6: CGTCGTCATTTCCGGCTTCAATTTC	28: TGCCCCGGCGTTTAGCGATCACAAC
7: GAGGAGAGCGAGTACTCGGGCTGCG	29: AAATACAGGCTCCACGACACGACCA
8: CGTGAAACCGCCCCAGCCTCGCCG	30: GGTTGCCCGCGCCCTTTTCCAGCC
9: ACTCGGAATCCCATGTGCTGACAGC	31: TCAGACAGGTTCCGCTCGATCAAGT
10: TCGACACCCGCTCTAGTTGACTTCC	32: GACCAAATAGGTATCGGCGTGTTC
11: GTGAGCAACGGCGCGGCAACCTGG	33: GACATGACGGCGGTGCCGCACTTGA
12: ATATCTGCTGCCCGCCCGGGAGAT	34: AAGTCACCTCGCCCACACCGTCGAA
13: GACCATCATTGCCATTCCCTCTCCC	35: TCCGTACGCTCGAAACGCTTCCAAC
14: GGTGTGATGCGGATGGTCGGCTCGG	36: CGAAATCCAGCACCACATCCGCAGC
15: CTTGAATAACGCGCAGTGAATTTCC	37: CGCGAACTCGTCCACAGTCCCCCTT
16: CGAGTTCCCGTCAGCGTCGTAAATC	38: CGTGGATGGCGGATGCGTTGTGCGC
17: GCGCCGGCCCGCGGGATGACTCCG	39: GACGATGGCCAGTAAATCGGCGTGG
18: CATGGACCCGGGCGAGCTGCAGATG	40: CGCCATCTGTGCCTCATACAGGTCC
19: TAACTGGCTTGGCGCTGATCCTGGT	41: GGAGCTTTCCGGCTTCTATCAGGTA
20: TTGACCTCGCCAGGAGAGAAGATCA	42: ATGGTGGGACATGGACGAGCGCGAC
21: TCGATGTCGATGTCCCAATCGTCGA	43: CGCAGAATCGCACCGGGTGC GGGAG
22: ACCGCAGACGGCACGATTGAGACAA	

¹⁾ The actual sequence in H37Rv of the underlined nucleotide is C.

²⁾ The actual sequence in H37Rv of the underlined nucleotides is GC.